



Complete Safety Solutions

SAVIOUR 'COVID 19 Kit' 'Basic'



General Description:

Protection of all body parts is a must specially while handling Fluids, aerosol particulates, pathogens or Live Virus , Bacteria etc. in any research Lab, Hospitals, Medical Centers, Nursing Homes etc. It is equally important when there is any infectious disease outbreak ,epidemic or Pandemic for complete protection of Medics, Para Medics and people at large. "Saviour" 'COVID-19 KITS' are just the right product for such applications as they are made up of PP Non woven Laminated System and doubly stitched to provide safety from ,micro-organisms , flying particulates ,direct or indirect and minute particles in the atmosphere or where spread is due to touch, contaminated blood, sputum or contamination in air etc.

MODEL NO: COVID 19 - BK

Specifications

The “Saviour” COVID-19 Kit **Basic** consists of the following;

- A) Full Body Coverall which is a Light weight disposable garment for protection against Virus made out of Non-Woven Laminated Fabric. The Coverall comes with a hood and has elastic at cuffs, waist & ankles having front zipper besides sleeves are provided with thumb pullers. The entire garment has sealed seams and comes in either white or sky blue colour.
- B) Disposable Hand gloves approved as per EN 374-2:2003 for Micro-organisms of minimum 2.5 to 3 MIL Thickness.
- C) 3 Ply Disposable Mask made out of PP non woven material.
- D) Clear Spectacles with wrap - around Polycarbonate lens. Single piece, ergonomically and sturdy design.
- E) Booties for extra protection of the feet.

Applications

- ◆ Control of Pandemics.
- ◆ Research Labs using Bacteria, virus.
- ◆ Disaster Management Authorities.
- ◆ Rescue and Relief Teams.
- ◆ Hospitals and Medical Centres.
- ◆ Pharmaceutical Industries.

Advantage

- ◆ Light weight, tough, puncture resistant Suit, Long Lasting with Extended life.
- ◆ Laminated hence protects against micro-organisms , flying particulates ,direct or indirect and minute particles in the atmosphere, Fluids, aerosol particulates, pathogens or Live Virus , Bacteria etc. Easily Disposable when contaminated.
- ◆ Hand Gloves for protection of Hands.
- ◆ Eyewear for Protection of eyes and super imposed with face Shield for protection of Eyes , Face and Mouth for extra Protection.
- ◆ 3 Ply Face Mask easily disposable.

Certification

- ◆ The Fabric Passes IS 1964:2001(Mass determination), ASTM E 96/ E 96M –05(MVT),BS EN 20811:1992(Resistance to Water Penetration), Passes ASTM F 1670 / F 1670M-08(2014)(Synthetic Blood Penetration Test). Also meets ASTM 2101-19 for Bacterial Filtration.
- ◆ Clear Spectacles approved as per EN-166:2001.
- ◆ Disposable Gloves meeting EN 374-2:2003 Standards.