



POWERED AIR PROTECTION SYSTEM [PAPS]

The 'PAPS' consists of a Complete Suit which is designed along with a hood provided with Optically Clear Wide Vision Visor and attached with Booties thus providing complete Body Protection. The System is supplied along with a Battery-Operated Powered Air Purifying system and Filter Cartridges.



COMPLETE PROTECTION

CERTIFICATION

P3 filters are used for protection against dust and toxic particulate according to The fabric is tested for Synthetic Blood Penetration Test Level 5 Exposure EN 143:1991, Class 3 standard and has a high filter efficiency (99.95%)

•Class 3 for Filter Cartridge.

•Fabric Tested to IS 7016 Part-I, II & IV.

•Fabric is tested as per ISO 2411, ISO 811

PANORAMIC VISOR
WITH 180° VIEW

THUMB LOOP FOR
BETTER & SECURE FIT
WITH ELASTIC WRIST

INTEGRATED BOOTIES

FLEXIBLE HOSE

EXHALATION VALVE

FILTER CARTRIDGE

POWERED AIR PURIFIER
SYSTEM (PAPS)



BENEFITS:



360° PROTECTION
& REUSABLE



LIGHT-WEIGHT
SOFT FABRIC



WATER
RESISTANT



SEAM SEALED

- Can be Autoclaved 100 % disinfection possible and can be safely reused.
- User Friendly as made from Soft Fabric hence comfortable to use.
- Provides complete protection as no other extra accessories to be used besides Hand Gloves and is also easy and quick to Don and Doff.
- Economical, Reusable and Made in India Product.
- The Whole System is light in weight and compact including the PA Purifier hence can be easily moved anywhere.