

Decontaminable Reusable Covid Suit

PVC FULL BODY COVERALL



General Description:

THE FULL BODY PROTECTIVE GARMENT 'HEAD-TO-TOE' (COVERALL) IS MADE OUT OF SPECIAL PVC MATERIAL WHICH IS TOTALLY WATERPROOF / RESISTANT TO FLUIDS.THE ENTIRE GARMENT IS TAPED AT SEAMS FOR COMPLETE PROTECTION AND IS ELASTICATED AT CUFFS (WITH THUMB LOOPS)AND ANKLES TO ENSURE TIGHT FITTING BESIDES HAVING ELASTICATED FACIAL FRONT IN HOOD.

THE FRONTAL OPENING IS PROVIDED FROM CHIN UPTILL BELOW THE CROTCH PORTION FOR EASE OF DONNING AND DOFFING AND CONSISTS OF A DOUBLE SLIDER ZIPPER WITH A COVERING FLAP FOR PROTECTION.