

SAVIOR KIOSK Covid-19 Test booth



MODEL NO:

KIOSK-19

Specifications

General Description:-

Saviour KIOSK Covid-19 Test Booth are for healthcare worker, These workers inside the kiosk collects the sample and then, follows the sanitisation process before proceeding to take the next sample. The collection process, fully contactless, gets over in five minutes. The main advantages of this procedure are that it needs fewer healthcare workers and strictly adheres to the norms of social distancing. Consisting of a small cubicle, it has PPE items such as gloves for the healthcare worker and a kit for collecting samples. It is also portable — it can be mounted on a vehicle and transported to any location. It can be particularly useful for collecting samples.

The kiosks have EN-374 approved built-in gloves, which can be used by health workers as protective shield while collecting samples of persons suspected with covid-19. All entry and exit to these kiosks operate using magnetic doors. They also have an exhaust fan, a chair and ultraviolet light for a patient's convenience. After each patient leaves the gloves are decontaminated with disinfectants and can be re-used.

Features:-

Aluminium Structure.

Dimensions: - Height 6 Feet X Weight 3 Feet X Length 4 Feet.

Front Dimensions 3 Feet x 3 Feet x 4 Feet Acrylic Window.

Side Infeed wall 2mm thickness rigid sheet along with the rubber bidding.

Floor Wooden Cover 12 mm with clean cover.