



The '**PAPS**' consists of a Complete Suit which is designed along with a hood provided with Optically Clear Wide Vision Visor and attached with Booties thus providing complete Body Protection. The System is supplied along with a Battery-Operated Powered Air Purifying system and Filter Cartridges.



CERTIFICATION

P3 filters are used for protection against dust and toxic particulate according to EN 143:1991, Class 3 standard and has a high filter efficiency (99.95%)

The fabric is tested for Synthetic Blood Penetration Test Level 5 Exposure

•Fabric Tested to IS 7016 Part-I, II & IV

• Fabric is tested as per ISO 2411 ISO 811

**PANORAMIC VISOR
WITH 180° VIEW**

FLEXIBLE HOSE

EXHALATION VALVE

**THUMB LOOP FOR
BETTER & SECURE FIT
WITH ELASTIC WRIST**

FILTER CARTRIDGE

INTEGRATED BOOTIES

**POWERED AIR PURIFIER
SYSTEM (PAPS)**

BENEFITS:



**360° PROTECTION
& REUSABLE**



**LIGHT-WEIGHT
SOFT FABRIC**



WATER RESISTANT



SEAM SEALED

- Can be Autoclaved 100 % disinfection possible and can be safely reused.
- User Friendly as made from Soft Fabric hence comfortable to use.
- Provides complete protection as no other extra accessories to be used besides Hand Gloves.
- Easy and quick to Don and Off.
- Economical, Reusable and Made in India Product.
- The Whole System is light in weight and compact including the PA Purifier hence can be easily moved anywhere.
- Battery Life: > 4 hours
- Weight: < 3 kg
- Spare Item : Cartridges (3 Nos.)

WARNING

- Check the product before use. Make sure it does not have any holes or is not damaged.
 - Do not use choked cartridges.
- Check the battery before use. Make sure it is fully charged.
- Once used, do not use PAPS again without decontaminating it.